

# MULTIPARA MONITOR F5



## Features

### Core

- ClassicARM9mainboardbasedonembeddedLinux
- Support storageof720hours trend tab leand graphre view,2 hours wave form review, 1000groups NIBP review and 200 alarm events review

### Body

- 12.1inchhighresolutionLCD
- Support display max. 11waveforms
- Support 7 channel ECG wave form display simultaneously

### Printer


- Built-in high-speed 50mm thermal printer Central System
- Wired or wire less connection
- Bed view

### Alarm

- Three-level acousto-optic alarm
- Human voice alarm
- Sensor-off alarm
- Paper out alarm
- Support alarm review
- Support alarm pause

### Linux OS

- Support data export for EMR connectivity by USB
- Multi-display mode
- NIBP self-test mode
- Support Data Export by USB
- 13 types of Arrhythmia analysis and real-time S-T segment analysis and pacemaker detection
- Drug calculation and titration table
- Hemo dynamic calculation
- Multi-language display




StandardDisplay



BigFontDisplay



ECGFullLead



TrendTable



NIBPTrend



OxyCRG

## Model Configuration

• **Standard Configuration** : 12.1inch LED, 5Lead ECG, NIBP, SpO2, Pulse Rate, Temperature, Respiration

• **Optional Configuration** : 2-Temperature, 2IBP, EtCO2, Masimo SpO2, Nellcor SpO2, Bed View

• **Optional Accessories** : Touch Screen, Printer, Central Monitoring Station, Neo/Ped Accessories, Wall-mounting, Trolley

**Optional** : EtCO2,Dual-IBP,12-leadsECG,Built-in thermal printer



**LIVO PRODUCTS INDIA (P) LTD.**

PLOT NO.20, FIRST FLOOR, BLOCK GI, LAWRENCE ROAD INDUSTRIAL AREA DELHI - 110035

PH: +91-11-43629200, 40196325

WWW.MEDICALPOINTINDIA.COM

EMAIL: MEDICAL\_POINT@YAHOO>COM