

Cardiotocograph Model: MCF-21K



Features

- Newly advanced A8 main board with Linux OS, more stable and safe
- Folded 10.4 inches high-brightness screen TFT LED, touch screen (optional)
- Supporting keyboard and mouse to input patient information
- 12-crystal pulse doppler transducer (waterproof for optional)
- FHR signal quality indicator helps optimize the probe position
- Large capacity storage space, 8640 hours memory for fetal traces
- Supporting window switch between standard interface and big font interface
- Event Marker makes it easy to annotate perceived fetal movements (FM) or other events on the strip chart

- Fetal Movement Profile (FMP) detects and automatically records gross fetal body movements
- Heart rate offset mode makes it easy to visually distinguish between the heart rates of twins by allowing you to offset the secondary FHR \pm 30BPM
- Cross Channel Verification (CCV) provide visual and audible indication when it automatically detects synchronous fetal or maternal heart rate/ pulse rate signals, indicating that you may be monitoring duplicate signals
- Fetal heart rate alarm provides both audible and visual alerts if the heart rate is outside of the user-defined high/low range.
- Easy-to-load Z-fold paper stacks itself neatly upon exiting the recorder, facilitating easy storage
- Internal rechargeable Li-on battery or AC power from external power supply.
- Build-in wireless network card, supporting wired or wireless connection to the central monitoring station
- Advanced Fischer NST Report (Fetal health report, optional)
- Support HL7 (Health Level Seven, optional)
- Optional Fischer No Stress Test (NST) report software allows interpretation of fetal heart rate and TOCO traces, and generates a printed report automatically for a reassuring NST

Fetal parameter					Maternal Parameter	Printer paper Width	Auxiliary Device			Advanced Fischer NST Report (Fetal Health Report)
Single FHR	Dual FHR	TOCO	AFM/MFM	Signal Overlap Alarm	NIBP, SpO ₂ , PR	112 mm	Lithium battery	Touch Screen/ Keyboard / Mouse/ Central Monitoring Station	Fetal Stimulator/ Waterproof Probe/ Wall-mount/ Trolley/ Bag	
√	√	√	√	√	x	√	√	Opt	Opt	Opt



UNICARE

SHENZHEN UNICARE ELECTRONIC CO., LTD.

3D, BLDG 4, TIAN AN INDL ZONE, NANYOU, NANSHAN, SHENZHEN, CHINA

Ph: +86-755-33369088

Website: www.unicare-med.com