

# MULTIPARA MONITOR XM750



**Suitable for Adult, Pediatric and Neonatal.**

**Reliable quality gives confidence, stable system makes perfect.**

## Product Feature

- Standard parameters: ECG, NIBP, RESP, PR, SpO<sub>2</sub>, TEMP
- Colourful and Clear 12.1" TFT screen, Backlight buttons
- Multiple display modes: Standard interface, Big Font, ECG standard full display, OXY Trend table, BP trend, View bed
- Ambulatory blood pressure technology, anti-movement
- Special design against high frequency surgical unit, and defibrillation protection
- Support Masimo SpO<sub>2</sub>-authorization partner 13 kind so arrhythmia analysis, 15 kinds of drug dose calculation
- Various languages
- Built-in detachable rechargeable lithium battery 4 hours battery life
- To connect CMS, other bed observation and software updating with wireless and wired mode
- Move able handle design Light weight portable, practical all-around, convenient operation

## Data Storage—stable and fast

- 8000 groups NIBP measurements
- 680 hours Trend data and Trend graphs
- 200 groups Alarm events reviewing
- 2 hours Wave forms reviewing
- Supporting for Arrhythmia analysis and drug dose calculation

## Alarms—more safe and reliable

- 3 level Audible and visual alarms

**Optional : EtCO<sub>2</sub>, Dual-IBP, 12-leads ECG, Touchscreen, Built-in thermal printer**



MEDICAL POINT



LIVO PRODUCTS INDIA (P) LTD.

PLOT NO.20, FIRST FLOOR, BLOCK GI, LAWRENCE ROAD INDUSTRIAL AREA DELHI - 110035

PH: +91-11-43629200, 40196325

WWW.MEDICALPOINTINDIA.COM

EMAIL: MEDICAL\_POINT@YAHOO.COM