

FETAL DOPPLER



FD-101CR



FD-201CR

Ultrasonic fetal doppler meets the fetus daily check and routine examination at home, clinic, community and hospital

Features

- Easy and convenient operation
- High-fidelity, crystal-clear sound
- Earphone and speaker are feasible
- High sensitivity doppler probe
- Low ultrasound dosage
- Display with white backlight LCD (**FD-101CR**) and color LCD (**FD-201CR**)

Technical Specification

- Ultrasound frequency: 2MHz
- Ultrasound intensity: <math><10\text{mW/cm}^2</math>
- Power supply: Alkalinity battery and Rechargeable battery
- Display: 35.04mm×28.03mm LCD
- FHR Measuring range: 50~230bpm
- FHR resolution: 1bpm
- FHR accuracy: $\pm 2\text{bpm}$
- Power consumption: <math><1\text{W}</math>
- Dimension :154mm ×103mm ×35mm
- Gross weight : approximately 500g

Configuration

- Main body
- Ni-MH battery (with adaptor)
- 2MHz probe

Functions	
Doppler type	Continuous
Built-in speaker	Yes
Battery indicator	Yes
Display	Digit
Various work mode	Yes
Battery	Ni-MH
Charger in device	Yes



LIVO PRODUCTS INDIA (P) LTD.

PLOT NO.20, FIRST FLOOR, BLOCK GI, LAWRENCE ROAD INDUSTRIAL AREA DELHI - 110035

PH: +91-11-43629200, 40196325

WWW.MEDICALPOINTINDIA.COM

EMAIL: MEDICAL_POINT@YAHOO>COM